AN

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Broad spectrum disinfectant

USE INSTRUCTIONS:

InSpec AN is designed for spraying, wiping and mopping applications. Hold spray approximately 15cm to 20cm from area to be treated. 5L versions, pour into an appropriate container for mopping. 100ml concentrate added to 2.4L of water is required for mopping. Apply to surface to ensure complete coverage.

TECHNICAL DATA:

Alcohol, chlorine and QUAT-free disinfectant cleaner. 4% RTU Solution pH 8.8-9.2, Concentrate pH 10.0.

STORAGE:

Store upright in original closed containers, away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is available in the Material Safety Data Sheet (MSDS).

SAFE HANDLING:

Always wear gloves and goggles or face protection. Always read the label and MSDS before use.

MATERIAL COMPATIBILITY:

Application of solutions, when used as directed, will not affect materials normally encountered in the cleanroom. See compatibility information in the Technical File.

BENEFITS:

- An exceptional non-oxidising biocide
- Bactericidal and Fungicidal
- Clean and disinfects
- 3 months validated 'in-use' shelf life
- Contents protected using 'bag-in-bottle' presentation with premium range.
- 1L and 5L formats manufactured with WFI quality water
- Multiple bags for cleanroom transfer
- Ideal rotational partner for InSpec QT, HA and OX
- Manufactured to GMP

TECHNICAL USER INFORMATION:

InSpec AN is a sterile disinfectant in pre-diluted and concentrated formats.

InSpec AN is manufactured to GMP in an ISO 6 cleanroom. The product is filled and capped in an ISO class 5 environment.

The solution is filtered through a 0.2 micron filter and gamma irradiated to a validated process (at 25-45kGy) to give a sterility assurance level SAL of 10⁻⁶.

Trigger sprays supplied with 900ml solution in a 1L container.

InSpec AN is supplied with Certificates of Conformity, Analysis, Irradiation and Sterility.

All formats of AN have a stabilised 24 month shelf life.

For screw cap bottles use the entire contents after opening.

Microbiological Minimum Efficacy Testing Contact

EN 1276 Bacteria	5 minutes	EN 13697 Bacteria	5 minutes
EN 1650 Fungi	15 minutes	EN 13697 Fungi	15 minutes

InSpec™ AN Product

Product	Pack Size	Case Size	Product Code
InSpec™ AN 900ml Trigger	1L trigger spray	6x 1L	ANWF120-900MLS
InSpec™ AN 5L Screw Cap	5L screw cap	2 x 5L	ANWF120-5LS
InSpec™ AN Concentrate	100ml screw cap	50 x 100ml	ANCNC20-100

Vietnam Official Distributor:

Decon Solutions Co., LTD

No 62 Road 64, Thanh My Loi Ward, District 2, HCM City, Vietnam
H: 09 0890 1955 **F:** +84 28 6287 3181 **E:** info@deconsolutions.com.vn
T: +84 28 6287 3180 - 38 143 585 **W:** www.deconsolutions.com.vn